

# Pilot Spray Dryers



### Features:

Aseptic GMP unit

Solvent / Aqueous Feeds

Co-Current Spray

Peristaltic Pump with forward / reverse and calibration facility

Twin cyclons (Optional)

Inert loop for Solvent Recovery (Optional)

High pressure jet cleaner (Optional)

### Technical specifications:

Models PSD-D / PSD-P	2 to 15 kg /hr water evaporation @ 250°C
----------------------	--

Inlet Air temperature	Ambient to 250°C
-----------------------	------------------

Heater	Electrical Heater
--------	-------------------

Feed Pump	Peristaltic Pump (Variable)
-----------	-----------------------------

Spray System	Two Fluid Spray Nozzle
--------------	------------------------

Power Supply	Three Phase, 440VAC, 50 Hz
--------------	----------------------------

\* Due to continuous improvement and development specifications / features / dimensions etc. may change without prior notice.

**Techno Search**  
**Instruments**  
technological essentials